

YZ30 Applanation Tonometer



The product is designed on the principle basis of Goldman tonometer. Two types are available: R and T type are suitable for Haag Streit type slit lamp.

Features:

- When used with slit lamp, it can be used to examine the eyes and measure the ocular pressure.
- Accurate measurement and total tolerance is less than 0.066KPa(0.55Hg).
- Directly get the ocular pressure and do not need to look up the conversion table.
- The measured ocular pressure is not affected by hardness of eye.
- The affected ocular volume is just 0.56 cubic millimeter.
- Adjustable measuring pressure ensures the long-term stability and reliability.

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

Web Site: www.66vision.com

E-mail: trade@66vision.com

Specifications:

Measuring Range	0~1064 Kpa
Light Ring Displacement	$1.53 \times 2 = 3.06\text{mm}$
Diameter of Prism Head	7mm
Moving Range of Prism Head	3mm

66VISION TECHCO.,LTD.

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

Web Site: www.66vision.comE-mail: trade@66vision.com